
# KerraCel® Ag

Gelling fiber silver dressing

## Instructions for use

### Description

KerraCel® Ag is a soft, sterile, nonwoven wound dressing made of sodium carboxymethylcellulose (CMC), cellulose fibers, and silver (0.2mg Ag/cm<sup>2</sup>) (1.7wt/wt%). The silver in the dressing provides an antibacterial barrier that inhibits bacterial growth in the dressing, as shown *in vitro*, for up to seven (7) days against Gram-positive and Gram-negative bacteria. KerraCel® Ag absorbs high amounts of wound fluid and creates a soft cohesive gel that intimately conforms to the wound surface and maintains a moist wound healing environment. KerraCel® Ag meets endotoxin limit specifications.

### Indications

Under the supervision of a health care professional, KerraCel® Ag may be used for the management of acute and chronic, partial and full thickness wounds including pressure ulcers, leg ulcers, diabetic foot ulcers, surgical wounds, traumatic wounds, first and second degree burns.

### Contraindications

Individuals with a known allergy or hypersensitivity to any of the components of KerraCel® Ag should not use the dressing.

### Precautions/Warnings

- **Caution:** Sterility is guaranteed unless pouch is damaged or opened prior to use. Do not resterilize
- KerraCel® Ag is not compatible with oil-based products, such as petrolatum
- Frequent or prolonged contact with the silver in the dressing may result in discoloration of the skin
- Avoid forceful removal and disruption of the healing wound
- If removing the dressing from the wound is difficult, the dressing should be fully saturated with sterile saline or water and removed slowly
- The dressing is not designed as a haemostat
- Dispose in normal waste if used in a household
- Dispose in clinical waste if used in a clinical establishment
- Store in a cool dry place (<77°F/25°C)
- The dressing may not be compatible with topical antimicrobials.
- Remove dressing prior to administering radiation therapy. A new dressing can be applied following treatment
- Avoid contact with electrode and conductive gels during electronic measurements (e.g. EEG and ECG)
- Silver is added to the dressing only as a preservative and to inhibit the growth of bacteria/microorganisms within the dressing

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Using KerraCel® Ag

1. Cleanse the wound in accordance with local clinical protocol.
2. Select an appropriate dressing size, so that the dressing overlaps the wound margin by approximately 1cm/0.5 inches.
3. When using KerraCel® Ag ribbon in tunneling wounds, leave at least 2.5cm/1 inch outside the wound for easy retrieval. Only pack deep wounds to 80% as KerraCel® Ag will expand to fill the wound space on contact with wound fluid.
4. For dry and minimally exuding wounds, place KerraCel® Ag on the wound and moisten with sterile saline or water over the wound area only.
5. The dressing should be secured with an appropriate moisture retentive secondary dressing (e.g. KerraFoam® Gentle Border, KerraLite Cool® Border, or other appropriate dressing).

### Frequency of Change

The dressing can be left in place for up to seven (7) days. However, this depends upon the condition of the wound and peri-wound skin. Check the dressing at regular intervals and change according to the clinical condition of the wound or when saturated. If the condition of the wound deteriorates, discontinue and contact a healthcare professional.

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KerraCel® is a registered trade mark of Crawford Woundcare Ltd.  
Ag Oxyalts™ is a trade mark of Exciton Technologies Inc.  
© Copyright Crawford Healthcare Ltd, 2016.

US Patent: 7300673, 7687076, 7998504  
Canadian Patent: 2,678,624

Manufactured for:  
**Crawford Healthcare Ltd**  
King Edward Court, King Edward Road, Knutsford, Cheshire, WA16 0BE, UK  
Tel +44 (0) 1565 654920 Email info@crawfordpharma.com [www.crawfordhealthcare.com](http://www.crawfordhealthcare.com)  
Distributed in the USA by:  
**Crawford Healthcare Inc.** 2005 S. Easton Road, Suite 203, Doylestown, PA 18901, USA  
Tel (855) 522-2211 [www.crawfordhealthcare.com/us](http://www.crawfordhealthcare.com/us)

KerraCel Ag Issue 01  
PN-04-0140  
Dec 2016  
Made in Canada