

Please keep these instructions, the black hard-plastic case, and the blue plastic return bag (with the affixed prepaid UPS shipping label) that was delivered with the product in a safe place.

Your supplier has partnered with The UPS Store for the return of KCI rental products. If you have any questions or need assistance with product returns, call your supplier or KCI at 1-800-275-4524.

When your therapy is complete:

**Step 1**

Call your supplier to discontinue rental of the therapy unit.

Remove the canister from the therapy unit. **DO NOT** return the canister to KCI; double-bag the canister and discard in the trash.



**Step 2**

Place *only* the therapy unit and power cords into the black hard-plastic case.

Discard the fabric carrying case, other items, literature, or documents containing your personal information, and any remaining unused or unopened disposables.



**Step 3**

Place the black hard-plastic case into the blue plastic return bag (with the affixed prepaid UPS shipping label).

Seal the blue plastic return bag.



**Step 4**

Drop the blue plastic return bag off at your nearest UPS Store; find a store near you at [www.ups.com/dropoff](http://www.ups.com/dropoff).

For a pickup call 1-800-PICK-UPS (1-800-742-5877); be prepared to provide UPS with the tracking information from the affixed prepaid UPS label on the blue plastic return bag.

