



ORDERING V.A.C.VIA™ NEGATIVE PRESSURE WOUND THERAPY

Ordering **V.A.C.Via™ Therapy**
is **streamlined** for you and
easy on your patients.



V.A.C.VIA™ THERAPY SYSTEM ORDERING INFORMATION

| Part Number | Description |
|----------------|---|
| VIKIT077D01/US | V.A.C.Via™ 7-Day Therapy System Kit |
| VIKIT07S05/US | V.A.C.Via™ Therapy Starter Kit 5-Pack |
| VIAGFS05 | V.A.C.Via™ GranuFoam™ Spiral Dressing, Small, 5-Pack |
| VIAGFM05 | V.A.C.Via™ GranuFoam™ Spiral Dressing, Medium, 5-Pack |
| VIACAN05 | V.A.C.Via™ Canister 250mL with Isolyzer 5-Pack |

For more information, visit www.KCI1.com or call **1-800-275-4524**

Device Disposal: At the end of therapy, follow local institutional protocols for infection control and waste disposal procedures for dressings and canisters. Local protocols for the disposal of the V.A.C.Via™ device should be based on applicable federal, state and/or local government environmental regulations for recycling electronic devices.



An Acelyty Company

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.

©2010 KCI Licensing, Inc. All rights reserved. All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and/or licensors. DSL#13-0653.US LIT29-B-236 (REV. 11/14)



V.A.C.VIA™ THERAPY SYSTEM PRODUCT SPECIFICATIONS*

| Feature | Description |
|---------------------------------|--|
| Unit Weight with empty canister | 0.7 lbs (0.32 kg) |
| Dimensions with empty canister | Length a. Length 6.4" (16.3cm); Width 3.5" (8.9cm); Depth 3.3" (7.0cm) |
| Pump | One diaphragm pump |
| User Interface | Membrane switch; 3-second press and hold |
| Therapy Delivery Modes | Continuous or Dynamic Pressure Control™ (DPC) Therapy |
| Therapy Pressure Options | Negative Pressure: -75mmHg or -125mmHg |
| Power Source | AC and battery |
| Battery Type | Custom lithium battery |
| Battery Life | ~ 9 hours |
| Battery Recharge Time | Less than 6 hours |
| Battery Indicator | Three-segment; low battery alarm when <2 hours |
| Canister | 250mL, omni-directional, with Isolyzer |
| Safety Alarms | System Fault, Blockage, Leakage, Low Battery, Therapy Life Exhausted |
| Electrical Data | External Power Supply Input: 100-240 VAC @ 50/60 Hz External Power Supply Output: nominally 9Vdc (18 Watt) Patient & Enclosure Leakage Current: < 100µA @ 132Vac @ 60 Hz |
| Electromagnetic Compatibility | Electromagnetic interference – Although this equipment conforms with the intent of the directive 89/336/EEC in relation to Electromagnetic Compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer. Radios, cell phones and similar devices may affect this equipment and should be kept at least 6.5 feet (2 meters) away from the equipment. Please use the guidelines and recommendations specified in Tables 204 and 206 of the EMC Insert. |
| IEC Classification | EC 60601-1 3rd edition Certified. *Type BF applied part *Class II equipment *IP22 ingress protection |
| Environmental Conditions | Transport and Storage Conditions: <ul style="list-style-type: none"> • Temperature Range: -20°F (-29°C) to 140°F (60°C) • Relative Humidity Range: 0-95%, non-condensing Operating Conditions: <ul style="list-style-type: none"> • Temperature: 41°F (5°C) to 104°F (40°C) • Relative Humidity: 15-93%, non-condensing • Altitude Range: Optimal Performance: -1,250ft to 9,850 ft. (-380m to 3,000m) |
| Packaging Dimensions | V.A.C.Via™ Therapy Starter Kit Shipper: 44 x 25.6 x 58.4 cm V.A.C.Via™ Therapy Starter Kit Each: 42.5 x 24.7 x 10.6 cm V.A.C.Via™ GranuFoam™ Spiral Dressing, Medium Shipper: 31.5 x 19 x 26.4 cm V.A.C.Via™ GranuFoam™ Spiral Dressing, Medium Pouch: 31.5 x 25.1 cm V.A.C.Via™ GranuFoam™ Spiral Dressing, Small Shipper: 31.5 x 19 x 26.4 cm V.A.C.Via™ GranuFoam™ Spiral Dressing, Small Pouch: 31.5 x 19 x 26.4 cm V.A.C.Via™ Therapy Canister Shipper: 46 x 20.8 x 14.6 cm |

For more
information, visit
www.kci1.com
or call
1-800-275-4524



*Specifications subject to change without notice